



SHIELD[®]

(Sterile Oral Surgery Kit with)



(ISO 9001:2015 & 13485:2013)

S.N.	Contents	Size	Qty
1	Oral Surgery Drape	80X80cm	1
2	Draw Sheet/Tray Cover	55X75cm	2
3	Chair Light Handle cover/Wrap	15X25cm	2
4	Surgical Handpiece & Suction hose sleeve		3
5	Patient Apron		1
6	Facemask		2
7	Surgeon Cap		2
8	Bouffant Cap		1
9	Surgical Instrument Tray (Plastic)		1
10	Surgical Gauze	05X05cm	4
11	High Suction Tip		1
12	Saliva Ejector		1
13	Rubber to hold I.V. line on surgical HP		4
14	Nylon ties to hold electric / I.V. cord		4
15	R.V.G. Sleeve		2
16	PVC Scrub Gloves		2 pair
17	Absorbent Towel	20X30cm	2



Technical Specification

This kit is made up of medical grade Non Woven, Hydrophobic Polypropylene SMS fabric as the main raw material. This material is preferred for its lightweight which is having 0.91 gravity only. It is lighter than that of 100% Hydrophobic PP fibers. It has excellent water repellency and air permeability resulting in comfort & 100% dryness.

Thus, provides resistance to Bacteria, Virus, Microbes, Fungi etc

Double medical grade packaging, **sterilized using Gamma Radiation** monitored, using protocol specified for surgical disposables, Lot & individual controls to ensure 100% sterility. (see color change of color indicator in individual

Instructions

1. Clean the unit and the working surface with 70% alcohol and disinfect the unit & working surface with disinfection spray.
2. Use sterile surgical kit for the patient & surface.

Recommendations

1. Mouth wash using chlorohexidine 2% for 30 seconds before starting the treatment.
2. Clean the unmasked portion of face by gauze soaked in chlorohexidine / Betadine.

Marketed by:

Sark Healthcare Pvt. Ltd.

J-8, 1st Floor, East Vinod Nagar Opp. Mayur Vihar-II, Delhi-110091

Email: info@sarkcare.com Website: www.sarkcare.com

Helpline: +91-9958007470