

# SINUS HELP KIT



Sark Healthcare



DOUBLE SAFE DEVICE

## SHK-01K SINUS HELP KIT

(Stopper : 7ea, Initial drill : 2ea, Diamond bur : 2ea, Hydraulic tip : 2ea, Mush : 1ea, lateral window bur : 2ea, window trephine : 1ea, window ring 1ea, Sinus curette 3ea.)

### COMPONENT

**STOPPER** : Diameter is just 4mm, smaller diameter than others.

**INITIAL DRILL** : Design for no bone burning and damage.

**DIAMOND DRILL** : Special diamond coated tip, 0.5mm tip stopper and shank stopper as well.

**HYDRAULIC TIP** : It can be fastened to osteotomy drill hole to maxillary sinus and lift sinus membrane without damage with water pressure.

**MUSHROOM** : Removing bone chipping and lifting sinus membrane safely.

**S WINDOW BUR** : It is designed for grinding bone for opening window safely with dual stopper and wide surface diamond.

**S WINDOW TREPHINE** : It is designed for making groove before using **S WINDOW RING** . ( 1mm blade )

**S WINDOW RING** : It is designed for grinding bone without damage of sinus membrane for open window.

It is better to use after **S WINDOW TREPHINE**

**SINUS CURETTE REGULAR** : 6 shape sinus curette considering the structure of human oral cavity and implant surgery.

### DETAIL VIEW



MU-2530  
MUSH

CRESTAL  
APPROACH

LATERAL  
WINDOW



SSR-40 SINUS SHORT REGULAR

SSR2-40 SINUS SIDE REGULAR

SLR-40 SINUS LONG REGULAR